

# Sepsis Risk Assessment Algorithm

## Complete Calculation Guide

### Purpose

This algorithm categorizes post-discharge patients into four risk tiers based on self-reported symptoms and home vital signs:

Tier	Label	Action
GREEN	Safe	Continue home monitoring
YELLOW	Concerning	Contact provider today
RED	Urgent	Seek immediate medical attention
RED_EMERGENCY	Emergency	Call 911

### Algorithm Flow (7 Steps)

INPUT: Patient survey responses

↓

STEP 1: Emergency Conditions → If triggered → RED\_EMERGENCY (exit)

↓

STEP 2: Hard RED Overrides → If triggered → RED\_EMERGENCY (exit)

↓

STEP 3: Detect Clinical Patterns → Determines pattern escalation level

↓

STEP 4: Calculate Base Score → Sum of individual symptom points

↓

STEP 5: Apply Context Bonuses → Additive adjustments from patient history

↓

STEP 6: Apply High-Risk Modifier → Score × 1.25 if applicable

↓

STEP 7: Determine Final Risk Level → MAX(score-based level, pattern-based level)

↓

OUTPUT: Risk tier + reasoning

### STEP 1: Emergency Conditions

**Logic:** If ANY of these are true, immediately return RED\_EMERGENCY. No further calculation.

Field	Condition	Result
fainted_or_very_dizzy	= true	RED_EMERGENCY

Field	Condition	Result
breathing_level	= 3 (extremely difficult)	RED_EMERGENCY
thinking_level	= 3 (not making sense)	RED_EMERGENCY
extreme_heat_or_chills	= true	RED_EMERGENCY
discolored_skin	= true (mottling/cyanosis/jaundice)	RED_EMERGENCY

## STEP 2: Hard RED Overrides

**Logic:** Critical vital sign thresholds that require immediate escalation. Evaluated after Step 1. No further calculation if triggered.

Field	Condition	Result
temperature_value	≥ 103.5°F	RED_EMERGENCY
oxygen_level_value	< 90%	RED_EMERGENCY
heart_rate_value	> 140 bpm	RED_EMERGENCY
breathing_level	= 3	RED_EMERGENCY
thinking_level	= 3	RED_EMERGENCY
blood_pressure_zone	= 3 (severe hypotension/crisis)	RED_EMERGENCY

*Note: breathing\_level = 3 and thinking\_level = 3 appear in both Step 1 and Step 2 as redundant safety nets.*

## STEP 3: Detect Clinical Patterns

**Logic:** Check for dangerous symptom combinations. Each pattern maps to an escalation level. Patterns can only **increase** risk, never decrease it. Priority order: RED\_EMERGENCY > RED > YELLOW.

### Infection Context Definition

Several patterns reference hasClearInfectionContext, which is true if ANY of:

- fever\_chills = true
- temperature\_zone ≥ 2
- urine\_appearance\_level ≥ 2
- wound\_state\_level = 3
- has\_cough = true AND mucus\_color\_level ≥ 2

### Pattern Definitions

**Pattern 1: Septic Shock** — only applies when blood\_pressure\_zone = 2; zone 3 is caught by the hard override in Step 2

```
IF blood_pressure_zone = 2 AND thinking_level ≥ 2 THEN:  
  IF hasClearInfectionContext = true → RED_EMERGENCY (exit immediately, skip  
  remaining patterns)  
  ELSE  
    → RED
```

### Pattern 2: Respiratory Failure

```
IF oxygen_level_zone ≥ 2 AND breathing_level ≥ 2 → RED
```

### Pattern 3: Silent Hypoxia

```
IF oxygen_level_zone ≥ 2 AND energy_level = 3 → RED
```

### Pattern 4: Septic Encephalopathy

```
IF fever_chills = true AND thinking_level ≥ 2 → RED
```

### Pattern 5: Compensated Shock

```
IF blood_pressure_zone ≥ 2  
  AND blood_pressure_systolic ≤ (baseline_bp_systolic OR 120)  
  AND heart_rate_zone ≥ 2  
→ RED
```

*This pattern only triggers for LOW blood pressure. A hypertensive crisis (>180) produces BP zone 3, which is caught by the hard override before this pattern is evaluated.*

### Pattern 6: Septic Hypothermia

```
IF temperature_value < 96.8°F → RED
```

### Pattern 7: SIRS Criteria

```
IF temperature_zone ≥ 2 AND heart_rate_zone ≥ 2 → YELLOW
```

### Pattern 8: Urosepsis

```
IF fever_chills = true AND urine_appearance_level ≥ 2 → YELLOW
```

Pattern 9: Cardiovascular Stress

```
IF energy_level = 3 AND heart_rate_zone ≥ 2 → YELLOW
```

**Output:** patternEscalation = highest level triggered (or null if none match)

STEP 4: Calculate Base Score

**Logic:** Sum points for each symptom present. Several fields have amplified scores when certain patient history flags are set.

Energy Level

Condition	Points
energy_level = 2 (fatigued)	+5
energy_level = 2 AND has_heart_failure = true	+10
energy_level = 3 (extremely fatigued)	+15
energy_level = 3 AND has_heart_failure = true	+25

Rationale: Fatigue in CHF patients signals cardiac decompensation earlier than in the general population.

Subjective Fever

Condition	Points
fever_chills = true AND has_thermometer = false	+10

Only scored when no thermometer is available. When a thermometer is present, temperature\_zone scoring applies instead.

Temperature

Condition	Points
temperature_zone = 2 (100–101.4°F)	+15
temperature_zone = 3 (<96.8°F or ≥101.5°F)	+40

Oxygen Level

Condition	Points
oxygen_level_zone = 2 (92–94%)	+15

Condition	Points
<code>oxygen_level_zone = 3 (&lt;92%)</code>	+40

### Subjective Heart Racing

Condition	Points
<code>heart_racing = true AND has_hr_monitor = false</code>	+5

*Only scored when no heart rate monitor is available.*

### Heart Rate (measured)

Condition	Points
<code>heart_rate_zone = 2 (101–120 bpm)</code>	+8
<code>heart_rate_zone = 2 AND has_heart_failure = true</code>	+15
<code>heart_rate_zone = 3 (&lt;60 or &gt;120 bpm)</code>	+30
<code>heart_rate_zone = 3 AND has_heart_failure = true</code>	+40

*Rationale: Tachycardia in CHF is a compensatory mechanism that signals decompensation earlier than in other patients.*

### Blood Pressure

Condition	Points
<code>blood_pressure_zone = 2 (20–39 mmHg below baseline)</code>	+8
<code>blood_pressure_zone = 3 (&lt;90, &gt;180, or 40+ below baseline)</code>	+20

### Mental Status

Condition	Points
<code>thinking_level = 2 (slow or foggy)</code>	+8

*`thinking_level = 3` triggers RED\_EMERGENCY in Steps 1 and 2 before scoring is reached.*

### Breathing

Condition	Points
<code>breathing_level = 2</code>	+8
<code>breathing_level = 2 AND (has_lung_condition = true OR has_heart_failure = true)</code>	+16

*Rationale: For COPD/asthma/lung fibrosis patients, mildly difficult breathing indicates exacerbation. For CHF patients, it may signal pulmonary oedema.*

## Urine Appearance

Condition	Points
<code>urine_appearance_level = 2</code> (cloudy/dark/smelly)	+15
<code>urine_appearance_level = 3</code> (very dark/bloody)	+60

## UTI Symptoms

Only scored when `has_recent_uti = true`.

Condition	Points
<code>uti_symptoms_worsening = 'same'</code>	+10
<code>uti_symptoms_worsening = 'same' AND has_urinary_catheter = true</code>	+18
<code>uti_symptoms_worsening = 'worsened'</code>	+30
<code>uti_symptoms_worsening = 'worsened' AND has_urinary_catheter = true</code>	+40

*Rationale: Catheter-associated UTIs (CAUTI) carry higher sepsis risk than uncatheterized UTIs.*

## Mucus / Cough

Condition	Points
<code>has_cough = true</code> (mucus level 1 or null)	+3
<code>mucus_color_level = 2</code> (yellow/green)	+10
<code>mucus_color_level = 2 AND has_lung_condition = true</code>	+18
<code>mucus_color_level = 2 AND has_recent_pneumonia = true</code>	+10 additional
<code>mucus_color_level = 3</code> (brown/pink/red)	+30
<code>mucus_color_level = 3 AND has_lung_condition = true</code>	+35
<code>mucus_color_level = 3 AND has_recent_pneumonia = true</code>	+10 additional

*Rationale: Colored or bloody mucus in patients with COPD, asthma, or lung fibrosis more reliably indicates exacerbation or superimposed infection. The pneumonia bonus stacks on top of the base or lung-condition score.*

## Wound Status

Condition	Points
<code>wound_state_level = 2</code> (looks different)	+10
<code>wound_state_level = 3</code> (infected: red/pus/swollen)	+60

## GI Symptoms

Condition	Points
<code>nausea_vomiting_diarrhea = true</code>	+5

---

## STEP 5: Apply Context Bonuses

**Logic:** Additive flat bonuses applied to the base score from Step 4. These reflect patient history risk factors independent of today's symptoms.

### Recent Discharge

Condition	Bonus
<code>days_since_last_discharge ≤ 7</code>	+15
<code>days_since_last_discharge 8–30</code>	+8
<code>days_since_last_discharge 31–90</code>	+3
<code>days_since_last_discharge &gt; 90</code>	+0

*Rationale: Sepsis recurrence risk is highest immediately post-discharge and decays over 90 days.*

### Prior Septic Shock

Condition	Bonus
<code>has_had_septic_shock = true</code>	+10

*Rationale: Prior septic shock is the strongest independent predictor of future severe sepsis. Even mild presentations warrant higher concern in this population.*

---

## STEP 6: Apply High-Risk Modifier

**Logic:** Certain patient groups receive a 25% score increase to compensate for atypical or attenuated symptom presentation.

### High-Risk Criteria

Patient is high-risk if ANY of these are true:

Field	Condition
<code>age</code>	<code>≥ 65</code>
<code>has_weakened_immune</code>	<code>= true</code>
<code>admitted_count</code>	<code>&gt; 1</code>
<code>on_immunosuppressants</code>	<code>= true</code>

*Rationale: Immunosuppressed patients may not mount a classic fever response; the score amplification compensates for attenuated symptoms. Patients admitted more than once have demonstrated recurrence risk.*

## Calculation

```
IF isHighRiskPatient = true:  
    totalScore = ROUND(baseScore × 1.25)  
ELSE:  
    totalScore = baseScore
```

---

## STEP 7: Determine Final Risk Level

**Logic:** Combine score-based thresholds with pattern-based escalation. Patterns can only increase the final level, never decrease it.

### Part A: Score-Based Level

```
IF totalScore ≥ 60:  
    scoreBasedLevel = RED  
  
ELSE IF totalScore ≥ 30:  
    scoreBasedLevel = YELLOW  
  
ELSE:  
    scoreBasedLevel = GREEN
```

### Part B: Pattern Escalation

```
IF patternEscalation = RED_EMERGENCY:  
    finalLevel = RED_EMERGENCY  
  
ELSE IF patternEscalation = RED AND scoreBasedLevel ∈ {YELLOW, GREEN}:  
    finalLevel = RED  
  
ELSE IF patternEscalation = YELLOW AND scoreBasedLevel = GREEN:  
    finalLevel = YELLOW  
  
ELSE:  
    finalLevel = scoreBasedLevel
```

## Final Formula

```
finalLevel = MAX(scoreBasedLevel, patternEscalation)
```



Where priority order is: GREEN < YELLOW < RED < RED\_EMERGENCY

## Zone Calculation Reference

### Temperature Zones

Zone	Range	Interpretation
1 (Green)	96.8°F – 99.9°F	Normal
2 (Yellow)	100°F – 101.4°F	Low-grade fever
3 (Red)	<96.8°F or ≥101.5°F	Hypothermia or high fever

Hard override: temperature ≥ 103.5°F → RED\_EMERGENCY (Steps 1–2, before scoring)

### Oxygen Zones

Zone	Range	Interpretation
1 (Green)	95% – 100%	Normal
2 (Yellow)	92% – 94%	Mild hypoxia
3 (Red)	<92%	Significant hypoxia

Hard override: oxygen < 90% → RED\_EMERGENCY (Steps 1–2, before scoring)

### Heart Rate Zones

Zone	Range	Interpretation
1 (Green)	60–100 bpm	Normal
2 (Yellow)	101–120 bpm	Tachycardia
3 (Red)	<60 or >120 bpm	Bradycardia or severe tachycardia

Hard override: heart rate > 140 bpm → RED\_EMERGENCY (Steps 1–2, before scoring)

### Blood Pressure Zones

Zone	Condition	Interpretation
1 (Green)	Within 20 mmHg of baseline	Normal
2 (Yellow)	20–39 mmHg below baseline	Moderate drop
3 (Red)	<90 mmHg, >180 mmHg, or ≥40 mmHg below baseline	Severe hypotension or crisis

Default baseline: 120 mmHg if *baseline\_bp\_systolic* not provided. Hard override: zone 3 → RED\_EMERGENCY.

# Complete Worked Examples

## Example A: Healthy Check-in

### Input:

- All emergency questions: false/normal
- Age: 45, no immune issues, admitted once
- Temperature: 98.6°F → zone 1
- Oxygen: 98% → zone 1
- Heart rate: 75 bpm → zone 1
- Blood pressure: 118 mmHg (baseline 120) → zone 1
- Energy: level 1, Thinking: level 1, Breathing: level 1
- All other symptoms: negative/normal

### Calculation:

1. Emergency conditions: none triggered
2. Hard RED overrides: none triggered
3. Patterns: none match
4. Base score: 0 (all vitals normal, no symptoms)
5. Context bonuses: none
6. High-risk modifier: no → totalScore = 0
7. Score-based level: GREEN (0 < 30)
8. Pattern escalation: none

### Final: GREEN

---

## Example B: SIRS Pattern (Fever + Tachycardia)

### Input:

- Temperature: 100.8°F → zone 2
- Heart rate: 112 bpm → zone 2
- All other vitals normal, no additional symptoms
- Age: 55, no high-risk flags

### Calculation:

1. Emergency conditions: none triggered
2. Hard RED overrides: none triggered
3. Patterns:
  - Pattern 7 (SIRS): `temperature_zone ≥ 2 AND heart_rate_zone ≥ 2` → YELLOW
4. Base score:
  - Temperature zone 2: +15
  - Heart rate zone 2: +8
  - **Base score = 23**
5. Context bonuses: none → adjusted score = 23
6. High-risk modifier: no → totalScore = 23

7. Score-based level: GREEN (23 < 30)
8. Pattern escalation: YELLOW → upgrades GREEN to YELLOW

**Final: YELLOW** — pattern escalated from GREEN

---

### Example C: Severe Hypotension (Hard Override)

#### Input:

- Blood pressure: 80 mmHg (baseline 120) → zone 3 (<90 threshold)
- All other vitals normal
- Thinking: level 1 (clear)

#### Calculation:

1. Emergency conditions: none triggered
2. Hard RED overrides: `blood_pressure_zone = 3` → **RED\_EMERGENCY** (exit)

**Final: RED\_EMERGENCY** — hard override triggered, no further steps

---

### Example D: Septic Shock Pattern

#### Input:

- Blood pressure: 95 mmHg (baseline 120) → zone 2 (25 mmHg below baseline)
- Thinking: level 2 (foggy)
- Temperature: 101.0°F → zone 2
- Heart rate: 105 bpm → zone 2
- No high-risk flags

#### Calculation:

1. Emergency conditions: none triggered
2. Hard RED overrides: none (BP zone = 2, not 3)
3. Patterns:
  - Pattern 1 (Septic Shock): `blood_pressure_zone = 2` AND `thinking_level ≥ 2`
    - Check `hasClearInfectionContext: temperature_zone ≥ 2` → **true**
    - → **RED\_EMERGENCY**, exit immediately

**Final: RED\_EMERGENCY** — septic shock pattern with confirmed infection context

---

### Example E: Multiple Symptoms, High-Risk Patient

#### Input:

- Age: 78 (high-risk: age ≥ 65)
- Temperature: 100.5°F → zone 2
- Heart rate: 108 bpm → zone 2
- Thinking: level 2 (foggy)

- Energy: level 2 (fatigued)
- Urine: level 2 (cloudy)
- Cough present, clear mucus (mucus level 1)
- GI symptoms present
- No `has_heart_failure`, no `has_lung_condition`

### Calculation:

1. Emergency conditions: none triggered
2. Hard RED overrides: none triggered
3. Patterns:
  - Pattern 7 (SIRS): `temperature_zone ≥ 2 AND heart_rate_zone ≥ 2` → YELLOW
4. Base score:
  - Temperature zone 2: +15
  - Heart rate zone 2: +8
  - Thinking level 2: +8
  - Energy level 2 (no CHF): +5
  - Urine level 2: +15
  - Cough, clear mucus: +3
  - GI symptoms: +5
  - **Base score = 59**
5. Context bonuses: none → adjusted score = 59
6. High-risk modifier: yes (age 78) →  $\text{totalScore} = \text{ROUND}(59 \times 1.25) = \text{ROUND}(73.75) = \mathbf{74}$
7. Score-based level: RED ( $74 \geq 60$ )
8. Pattern escalation: YELLOW (RED > YELLOW, no upgrade)

**Final: RED** — score-based level wins over pattern

---

### Example F: Compensated Shock (Tachycardia + Hypotension)

#### Input:

- Blood pressure: 98 mmHg (baseline 125) → zone 2 (27 mmHg below baseline)
- Heart rate: 118 bpm → zone 2
- Thinking: level 1 (clear)
- No fever, no infection signs, no high-risk flags

### Calculation:

1. Emergency conditions: none triggered
2. Hard RED overrides: none triggered
3. Patterns:
  - Pattern 5 (Compensated Shock): `blood_pressure_zone ≥ 2 AND bp_systolic (98) ≤ baseline (125) AND heart_rate_zone ≥ 2` → RED
4. Base score:
  - BP zone 2: +8
  - Heart rate zone 2: +8
  - **Base score = 16**

- 5. Context bonuses: none → adjusted score = 16
- 6. High-risk modifier: no → totalScore = 16
- 7. Score-based level: GREEN (16 < 30)
- 8. Pattern escalation: RED → upgrades GREEN to RED

**Final: RED** — pattern escalated from GREEN

*Clinical rationale: The heart is racing to compensate for falling blood pressure — a pre-shock state requiring urgent evaluation even though individual scores are low.*

### Summary Table: All Point Values

Category	Condition	Points
Energy	Level 2 (fatigued)	+5
Energy	Level 2 + <code>has_heart_failure</code>	+10
Energy	Level 3 (extremely fatigued)	+15
Energy	Level 3 + <code>has_heart_failure</code>	+25
Subjective fever	<code>fever_chills</code> + no thermometer	+10
Temperature	Zone 2	+15
Temperature	Zone 3	+40
Oxygen	Zone 2	+15
Oxygen	Zone 3	+40
Subjective HR	<code>heart_racing</code> + no monitor	+5
Heart rate	Zone 2	+8
Heart rate	Zone 2 + <code>has_heart_failure</code>	+15
Heart rate	Zone 3	+30
Heart rate	Zone 3 + <code>has_heart_failure</code>	+40
Blood pressure	Zone 2	+8
Blood pressure	Zone 3	+20
Thinking	Level 2	+8
Breathing	Level 2	+8
Breathing	Level 2 + <code>has_lung_condition</code> or <code>has_heart_failure</code>	+16
Urine appearance	Level 2	+15
Urine appearance	Level 3	+60
UTI symptoms	Unchanged ( <code>has_recent_uti</code> )	+10

Category	Condition	Points
UTI symptoms	Unchanged + <code>has_urinary_catheter</code>	+18
UTI symptoms	Worsened ( <code>has_recent_uti</code> )	+30
UTI symptoms	Worsened + <code>has_urinary_catheter</code>	+40
Cough	Present, clear/no mucus	+3
Mucus	Level 2 (yellow/green)	+10
Mucus	Level 2 + <code>has_lung_condition</code>	+18
Mucus	Level 2 + <code>has_recent_pneumonia</code>	+10 additional
Mucus	Level 3 (brown/pink/red)	+30
Mucus	Level 3 + <code>has_lung_condition</code>	+35
Mucus	Level 3 + <code>has_recent_pneumonia</code>	+10 additional
Wound	Level 2 (different)	+10
Wound	Level 3 (infected)	+60
GI symptoms	Present	+5
<b>Context bonus</b>	Discharged $\leq 7$ days ago	+15
<b>Context bonus</b>	Discharged 8–30 days ago	+8
<b>Context bonus</b>	Discharged 31–90 days ago	+3
<b>Context bonus</b>	Prior septic shock history	+10

## Summary Table: Risk Thresholds

Condition	Result
Emergency question triggered (Step 1)	RED_EMERGENCY
Hard vital sign override triggered (Step 2)	RED_EMERGENCY
Pattern = RED_EMERGENCY	RED_EMERGENCY
<code>totalScore <math>\geq</math> 60</code>	RED
Pattern = RED	RED
<code>totalScore <math>\geq</math> 30</code>	YELLOW
Pattern = YELLOW	YELLOW
Otherwise	GREEN